

NATIONAL UNIVERSITY OF SINGAPORE

Doc no: CIBA/RA/Eq/015 Experiment-Based Risk Assessment Form

Name of Department PHYSICS Location of Lab S11-02-09
 Name of Laboratory CIBA Optical Materials & Devices Lab Name of PI Asst Prof Andrew Bettiol
 Name of Researcher/LO Yan Yuanjun; Yang Chengyuan; Vanga Sudheer Kumar; Shuvan Prashant Turaga Name of Activity/Experiment Two Photon Lithography

No	Description/Details of Steps in Activity	Hazards	Possible Accident / Ill Health & Persons-at-Risk	Existing Risk Control (Mitigation)	Severity	Likelihood (Probability)	Risk Level	Additional Risk Control	Person Responsible	By (Date)
1	Laser alignment	1. Laser light(Class 4)	1. Laser light shining directly into eyes can cause permanent blindness. (burning of corneas)	1. All users should wear laser goggles of appropriate wavelength. 2. Black non-reflective boards barricading areas where laser beams are aligned. 3. Keep the MIRA cavity closed at all times. 4. Use IR card for alignment.	3	1	3			
2		2. Fire hazard	1. Fire caused by high power femtosecond infrared light 2. Skin burn due to laser light	1. Always close the laser shutter when the laser is not in use. 2. Reduce the power to minimum during laser alignment. 3. Cover the light path when laser is in use. 4. No paper or anything that can catch fire in the light path 5. Use metal shields to block relected light	2	1	2			
3		3. Reflected laser light	1. Reflected laser light can cause permanent blindness.	1. No jewellery allowed when working with lasers. 2. 'LASER IN USE' sign lighted when laser work is carried out.	1	1	1			
		4. Overloading of laser power source	1. Electrical short circuit	1. Label stating maximum allowed current allowed for each power supply. 2.Do not open power supply	1	1	1			
	Using Two Photon Optical Setup	Fire hazard	Focused beam can burn paper or flammable materials	1. Don't use paper for aligning the beam 2. Use IR Card for alignment 3. Close the laser shutter when not using the beam for fabrication.	2	1	2			
		Reflection Hazard	The reflection from the sample or a glass slide placed underneath objective	1. Don't align the beam at eye level. 2. Wear appropriate Optical Goggles while doing the alignment. 3. Keep the beam power low for alignment.	1	1	2			

Conducted By Yan Yuanjun
Yang Chengyuan
Vanga Sudheer Kumar
Shuvan Prashant Turaga

Approved By Name Asst Prof Andrew Bettiol
 Signature _____
 Approval date 1-Nov-11 Next Revision date 1-Nov-14
 (Maximum 3 years)