
 <b>NUS</b> National University of Singapore  <b>CIBA</b> Faculty of Science, Dept of Physic , Centre of Ion Beam Applications	Procedure No:	CIBA/SOP/Eq /003
Title:  <b>Spin Coater</b>	Rev No:	0002
	Issue Date:	31 Oct 2011
Prepared by: Liu Fan & Wang Yinghui	Page:	Page 1 of 3
	Approved by: Asst Prof Jeroen Van Kan	Review Date: 04 Nov 2011

## 1 Objective:

This Standard Operation Procedure states how Spin Coater should be operated.

## 2 Responsibilities:

### 2.1 Director / HOD / PI

The Director/HOD/PI has overall responsibility for ensuring a system is established for the safe use of the spin coater.

### 2.2 Designated Person

There shall be a designated person to oversee the correct operation and maintenance of the spin coater.

a. This person shall periodically inspect the spin coated to ensure its operational performance.

b. He/she will make necessary arrangements for repair works of the spin coater.

d. He/she will report to the Director/HOD/PI unsafe practices by spin coater users.

### 2.3 Staff/ Research personnel

a. Spin coater users shall attend appropriate training on the safe use of the machine.

b. Users shall report any injuries, defects or breakdowns to their supervisor.

## 3 Procedures:

3.1 Switch on power supply and Transformer (black box)

3.2 Switch on Vacuum Pump

3.3 Turn Lever Towards the Gold colored bolt (for the N<sub>2</sub> gas Make sure N<sub>2</sub> ≥ 60 PSI)

3.4 Use aluminum foil to cover up the inside of the coater

3.5 Set Program


## 4 Programming Procedures (PGM mode):

4.1 . Press **Program Select** button to select program (program A to T  
\*\* Don't use C this is for cleaning!)

4.2 . Press **F1** button to switch between "program" mode and "off" mode (indicated as PGM and OFF on the interface respectively)

4.3 . In PGM mode, press **STEP** button to select the step to be programmed

4.4 . Press ← or → buttons to switch between the set data (RPM, time etc) to be changed

 <b>NUS</b> National University of Singapore <b>CIBA</b> Faculty of Science, Dept of Physic , Centre of Ion Beam Applications		Procedure No:	CIBA/SOP/Eq /003
Title:  <b>Spin Coater</b>		Rev No:	0002
		Issue Date:	31 Oct 2011
		Page:	Page 2 of 3
Prepared by: Liu Fan & Wang Yinghui	Approved by: Asst Prof Jeroen Van Kan	Review Date: 04 Nov 2011	

4. 5 . Press ↑ or ↓ buttons to increase or decrease the values of the set data
4. 6 . To delete steps, press **DEL STEP** button and press **ENTER** button
4. 7 . ACL (on the interface) is the acceleration index value and be adjusted while the number to the right of it is the calculated acceleration (RPM/second) which is not adjustable but is dependent on the selected ACL

## 5.1 Running the Program:

5. 1 Load sample onto chuck  
**Make sure the Sample is larger than the chuck**
5. 2 Put resist      **\*\*MAKE SURE YOU DON'T PUT ANY CHEMICALS INTO THE CHUCK!**
5. 3 Close lid and Press **VACUUM** button
5. 4 Press F1 button and ensure that coater is in "OFF" mode
5. 5 Press **RUN/STOP** button and program will start running (Coater will be in "RUN mode as indicated on the interface)
5. 6 When program ends, press **VACUUM** button, open lid and take out sample

## 6 Shut down Procedures

6. 1 Clear aluminum foil
6. 2 Cover chuck with Petri Dish and drop few drops of acetone on the dish
6. 3 Press **VACUUM** button and close the lid and run program "C"
- 6. 4 Remove any resist on the chuck with acetone.**
6. 5 When program ends, open lid and wipe inside of Coater with acetone and dry it taking care not to let any chemicals enter the vacuum path via the chuck
6. 6 Press **VACUUM** button
6. 7 Store the Petri Dish
- 6. 8 Switch off vacuum pump**
- 6. 9 Switch off Transformer and power supply**
6. 1 0 Wait 2 mins after last use before switching off the N2.

Title:

## Spin Coater

Rev No:

0002

Issue Date:

31 Oct 2011

Page:


Page 3 of 3

Prepared by:  
Liu Fan & Wang Yinghui

Approved by:  
Asst Prof Jeroen Van Kan

Review Date:  
04 Nov 2011

### 7 Summary guide:



### 8 References:

Manuals kept on the shelves at S7-01-01A near grey area cabinet.

### 9 Appendix:

Nil